

New !!

New !!

New !!

New !!

EASY CAP EEG-Recording Caps

are well known for

- Superior Signal Quality from Ag/AgCl Sintered Electrodes
- Lowest Impedances together with highest Subject Comfort
- Customer-specific Electrode Arrangements

Serving prime research labs now for 8 years we compiled a

Ready-to-use EEG Recording Cap Set for 64 / 80 Channels

**incorporating all 10%-system electrode positions
for all EEG machines with 64 or 80 Channels.**

This package contains everything you need
to start recording right out of the box:

- Cap with 74 equidistant electrode positions (10%-system)
- 78 Ag/AgCl sintered electrodes, number labeled 1 – 78
- 2 Ear-lobe electrodes
- Chest belt set
- All necessary consumables

and as a **totally new feature** to shorten mounting time

- **a quick-connect interface!**

To get all this just order

EC80-#

(Replace the # by the preferred cap size)



EASY CAP EC80

Standardized Ready-to-use EEG Recording Cap Set for 64/80 Channels

consisting of

- | | | |
|----|-----------|--|
| 1 | SU-# | Subtemporal cap for # cm head circumference*
<i>equipped with 74 electrode adaptors A06 according to Montage No. 1 (see graphic)</i> |
| 78 | B10-S-120 | AgCl sintered ring electrode
<i>with 120 cm Heavy-duty leadwire and 1,5mm touchproof safety socket</i> |
| 1 | C4 | Pair of Nihon-Kohden ear clips
<i>with appropriate electrodes</i> |
| 2 | KV | 38-Channel-Electrode board adaptor |
| 2 | KI | 38-Channel-Input box |
| 1 | GS | Chest belt set:
<i>1 Chest belt, 2 chin straps, 1 chin belt</i> |
| 1 | V17 | Abralyt 2000 (1 kg)
<i>Chloride-free, abrasive electrolyte gel</i> |

Consumables:

Alcohol, cotton swabs, syringe, measuring tape, spare adaptors, cleaning brush, manual, channel layout graphic

The catalogue number of the entire package is **EC80-#***.

When re-ordering, please use the catalogue numbers of the single items, e.g.


- | | | |
|---|------|--|
| 1 | SU-# | Subtemporal cap for # cm head circumference*
<i>equipped with 74 electrode adaptors A06 according to Montage No. 1</i> |
|---|------|--|

Our full product list can be seen at www.advancedmedicalequipment.com

* Cap sizes are given in cm head circumference:

- Adult caps: 62, 60, 58, 56, 54
- Children caps: 54, 52, 50
- Infant caps: 36, 38, 40, 42, 44, 46, 48

Montage No. 1
74 Equidistant electrode sites, representing full 10%-system
Spare electrodes can be placed on skin for use as EOG, ECG, etc.



Presented by:

advanced
MEDICAL
equipment ltd

Tel: 01403 260156

Fax: 01403 260175

info@advancedmedicalequipment.com