


AUDIOMETRY ACCESSORIES



■ Insert Earphones for Audiometry

Kit Includes: Earphones 10 Ohm, Dual-mono 7' cable assembly, 50 foam eartips (regular, 13 mm), 50 foam eartips (baby, 10 mm), 2 foam eartips (jumbo, 18 mm) and Velcro

Product Description

STIM or STIM2 10 Ohm

insert earphones 00100268

Part Number



■ Disposable Foam Ear Inserts

Foam eartips developed for ER-3A insert earphones are produced with constant dimensions to ensure proper calibration and test accuracy. To accommodate differences in ear canal size, there are three sizes of foam eartip (small, medium and large).

Product Description

10mm foam ear inserts 50/pkg. 97000127

13mm foam ear inserts 50/pkg. 97000128

18mm foam ear inserts 24/pkg. 97000129

Part Number



■ Replacement Front Tubes with Adapters

Product Description

Includes 4 tubes and 4 adaptors 97000130

Part Number



■ Sound Tube Adapters

Product Description

Package of 10 front tube adaptors ... 97000131

Part Number