

DISINFECTANTS AND CLEANERS



■ Envirocide Disinfectant

Envirocide is a convenient, ready-to-use, intermediate level disinfectant which is effective against TB, HBV, HCV, viruses (hydrophilic and lipophilic), bacteria (including MRSA and VRE) and fungi.

Product Description	Part Number
Envirocide 24 oz. spray bottle	00104209
Envirocide 1 gallon bottle	00104210



■ Control III Disinfectant

Formulated for use after dilution on pre-cleaned surfaces and equipment. Non-irritating on skin, and has a pleasant odor. A 10 minute contact time is recommended for disinfection.

Product Description	Part Number
Control III Gallon (ready to use)	00103724
Control III Concentrate Pint	00103725
Test Strips 15/Pack	00103726



■ Control III Elite

A ready to use EPA approved detergent and disinfectant in one. Elite is effective against TB as well as other bacteria and fungi.

Product Description	Part Number
Control III Elite Gallon	00103727
Control III Elite 32 oz. Spray bottle	00103728



■ CaviWipes

Presaturated with CaviCide, CaviWipes are proven to kill both TB (in 5 minutes) and HBV (in 3 minutes) with quick, easy-to-use, timesaving convenience. CaviWipes are durable, nonwoven, nonabrasive towels that won't bunch up during use. They offer better surface contact than other wipes, requiring less force to apply solution to a surface.

Product Description	Part Number
Caviwipes XL 12" x 10" 66 count can ...	00103731
Caviwipes 6" x 6" 160 count can	00103732



■ CaviCide Disinfectant

CaviCide is a convenient, ready-to-use, intermediate level disinfectant.

Product Description	Part Number
CaviCide 24 oz. spray bottle	00103729
CaviCide 1 Gallon bottle	00103730



■ T-Spray Disinfectant

Cleans and disinfect ultrasound transducers, probes, sensors, and other hard surfaces. Simply spray covering the surface with T-Spray and wipe with a clean cloth.

Product Description	Part Number
T-Spray 60ml bottle	00103443
T-Spray 250ml bottle	00103444



■ Dial Antibacterial Hand Soap

Remove as much as 97% of contaminating bacteria and viruses with Liquid Dial Soap. Dial is mild to the skin, even with frequent washing, and it provides broad antimicrobial effectiveness.

Product Description	Part Number
Dial 7.5 oz. pump bottle	95000061
Dial 1 Liter refill cartridge	95000062