

ELECTRODE APPLICATION



Adhesive Disks

Two sided adhesive disk holds electrodes in place.
Sold 100 per pack.

Product Description	Part Number
4mm ID x 19mm OD	97000091
8mm ID x 22mm OD	97000093
12mm ID x 32mm OD	97000092



V-Shaped Electrode Washer

Single-sided adhesive washer secure electrodes and sensors. Sold 100 per pack.

Product Description	Part Number
V-Shaped adhesive washer	97000090



Adhesive Foam Disks

Foam disks with single-side adhesive.
Sold 100 per pack.

Product Description	Part Number
Adhesive foam disks	97000089



Cotton Tip Applicator

6" non-sterile wooden cotton tip applicator.
Sold 1000 per pack.

Product Description	Part Number
6" non-sterile applicator	95000009



Alcohol Prep Pads

Pads measure 2 3/8" x 1 1/8" and come individually packed. 200 pads per pack.

Product Description	Part Number
200 pack alcohol prep pads	95000013
16 oz. bottle 70% Isopropyl alcohol	95000012



Alcohol Prep Pads with Pumice

Cleanse skin in preparation for electrode placement. Pumice removes excess oil for better electrode to skin contact. Individually packaged 100 per box.

Product Description	Part Number
Prep pad with Pumice	95000016



Adhesive Remover Pad

Effectively softens and removes tape and residue from skin. 100 pads per box.

Product Description	Part Number
Adhesive remover pads	95000011



■ Cotton Balls

Medium size cotton balls. 200 per bag.

Product Description	Part Number
Cotton Balls	95000024



■ Gauze

Non-sterile multi-purpose cotton gauze squares. 200 per pack.

Product Description	Part Number
2" x 2" 200 pack	95000025
4" x 4" 200 pack	95000026



■ Coban Self-Adherent Wrap

Elastic wrap is versatile and easy to apply. Easily adheres to itself without the need for any pins, clips, or tape.

Product Description	Part Number
2" x 5 yds. 36 per box	95000022
1" x 5 yds. 20 per box	95000023



■ Head Net Tubular Bandage

A stretchable tubular elastic net material used to help hold electrodes in place. Comes in a 25 yard box and is easily cut to length.

Product Description	Part Number
Size 5 tubular net – Small adult/child ...	95000046
Size 6 tubular net – Adult/Large	95000048



■ China Marker

Water soluble red china marker. 12 per box.

Product Description	Part Number
Red China Marker	00101840



■ Tape Measure

Self retracting easy to read tape measure. Both metric and inch scales.

Product Description	Part Number
Tape measure	95000038



■ Cover Roll Stretch Tape

Non-woven stretch adhesive bandage is excellent for areas requiring flexibility while providing superior holding power. Hypoallergenic, air permeable, and latex-free. Ideal for sensitive skin. Sold by the roll.

Product Description	Part Number
2" x 10 yds. Roll	00103467



■ Micropore Surgical Tape

Microporous paper tape is comfortable and lightweight. Holds even when wet. Hypoallergenic and latex-free.

Product Description	Part Number
1/2" x 10 yds. 24 per box	95000035
1" x 10 yds. 12 per box	95000032
2" x 10 yds. 6 per box	95000060



■ Transpore Surgical Tape

Perforated for easy tearing and porous for greater breathability. Hypoallergenic.

Product Description	Part Number
1/2" x 10 yds. 24 per box	95000034
1" x 10 yds. 12 per box	95000031



■ Durapore Surgical Tape

Strong, highly adhesive hypoallergenic silk tape features maximum comfort and easy tearing for quick application.

Product Description	Part Number
1" x 10 yds. 12 per box	95000030
2" x 10 yds. 6 per box	95000063



■ Blenderm Surgical Tape

Transparent and waterproof. Conforms to body areas which are typically hard to tape. Hypoallergenic clear plastic minimizes the risk of skin irritation. Latex-free.

Product Description	Part Number
1" x 5 yds. 12 per box	95000033



■ Latex exam gloves

Lightly powdered gloves offer a good grip. Beaded cuffs offer greater protection. Non-sterile, 100 per box.

Product Description	Part Number
Small latex exam gloves	95000027
Medium latex exam gloves	95000028
Large latex exam gloves	95000029



■ Latex-free Purple Nitrile Exam Gloves

Nitrile gloves provide extra coverage to protect infectious liquids. Their purple color helps users immediately distinguish them from other latex or synthetic options. Non-sterile, powder-free, latex-free, extra long cuffs. 50 per box

Product Description	Part Number
Small powder-free/Latex-free exam gloves	00103602
Medium powder-free/Latex-free exam gloves	00103603
Large powder-free/Latex-free exam gloves	00103604