

EMG AND IOM NEEDLES AND ACCESSORIES

■ Neuroline® Disposable Monopolar EMG Needle



Neuroline Monopolar needle electrodes are made of stainless steel and feature the optimum combination of strength and flexibility. Each needle is equipped with a color coded ergonomic hub which provides ease of identification and better control. Needles are available with pre-attached 30" lead wire or with no lead wire. Needles without attached lead wire require reusable adaptor cable. Needles are sold 40 to a box.

Product Description	Part Number
25mm (1") x .36mm (28G) Needle with Pre-attached Lead wire	97010003
38mm (1.5") x .36mm (28G) Needle with Pre-attached Lead wire	97010004
38mm (1.5") x .45mm (26G) Needle with Pre-attached Lead wire	97010022
50mm (2") x .45mm (26G) Needle with Pre-attached Lead wire	97010005
75mm (3") x .50mm (25G) Needle with Pre-attached Lead wire	97010023

Product Description	Part Number
Reusable Monopolar cable	97010010

■ Neuroline® Inoject™ Disposable Botox Needle



Single patient EMG Needle electrodes for injections with Botulinum Toxin. Available in six different sizes to meet your needs. Each size is manufactured with the optimum combination of strength and flexibility. Each needle is made of stainless steel with a pre-attached 30" Leadwire. The needle hub has a Luer Lok fitting, which secures the needle to the syringe. Needles are sold 10 to a box.

Product Description	Part Number
25mm (1") x .30mm (30G) Needle	97010011
30mm (1.2") x .36mm (28G) Needle	97010012
35mm (1.4") x .40mm (27G) Needle	97010013
38mm (1.5") x .45mm (26G) Needle	97010014
50mm (2") x .50mm (25G) Needle	97010015
75mm (3") x .55mm (24G) Needle	97010016