

ELECTRODES - EEG, EOG, ECG

■ Genuine Grass® Brand Electrodes

Get unsurpassed recording clarity, long life and dependability right from the start. Every Grass® electrode is made of the finest materials and is precision handcrafted in an exclusive 12-step manufacturing process. The result is superior quality that assures reliable, accurate EEG, EMG, ECG and EOG recordings.

Each package contains 10 electrodes.



Gray Connectors

Product Description	Part Number
Gold Cup	7000-0042-00
60" 10mm Gold Cup	7000-0043-00
72" 10mm Gold Cup	7000-0044-00
96" 10mm Gold Cup	00103445
120" 10mm Gold Cup	00103446

Genuine Grass Electrodes with No-Tangle Wires

Product Description	Part Number
48" 10mm Gold Cup	00103452
60" 10mm Gold Cup	00103453
72" 10mm Gold Cup	00103454
96" 10mm Gold Cup	00103455
120" 10mm Gold Cup	00103456

■ Reusable 10mm EEG Cup Electrodes

Available in Gold, these electrodes are manufactured with a superior molded connector. These reusable electrodes are Teflon coated tinsel wire providing a stronger and lighter wire when compared to conventional electrodes. Each package contains 10 –10mm cup electrodes that terminate into the standard 1.5mm connector.



Product Description	Part Number
48" 10mm Gold Cup	00103459
60" 10mm Gold Cup	97000050
72" 10mm Gold Cup	00103460
96" 10mm Gold Cup	00103461

■ Reusable Cast Disc Electrodes

Reusable Gold/Silver Silver Chloride (Ag/AgCl) 10mm diameter disc electrodes with Polyurethane "tangle free" lead wire and molded 1.5mm connector.

Each Package contains 10 electrodes.

Product Description	Part Number
1.5m 10mm Gold Disc	97000057
1.2m 10mm Ag/AgCl Disc	97000055
1.5m 10mm Ag/AgCl Disc	97000059



■ Reusable Silver/Silver Chloride Electrodes

Each electrode is made up of a homogenous mixture of silver-silver chloride compressed and sintered to strengthen and stabilize. These electrodes never require rechloriding and are best applied to the skin using adhesive washers.

Sold individually

Product Description	Part Number
1.2m Sintered silver/silver chloride electrode	97000152
1.5m Sintered silver/silver chloride electrode	97000153



■ Disposable Cup Electrodes

Single patient electrodes feature high quality silver/silver chloride sensor and tangle free lead wires. These electrodes make it possible to use new, unused electrodes with each patient which reduces the risk of cross-contamination.

Product Description	Part Number
100cm (40") Pack of disposable electrodes	00103436
150cm (60") Pack of disposable electrodes	00103437





Reusable Earclip Electrodes

This product consists of two EEG electrodes mounted as an ear clip. Electrodes are flat for easy application to the earlobe. Product is sold in pairs.
Earclip Electrode with 6cm leadwire

Earclip Electrode with 1.5m (60") Leadwire

Product Description	Part Number
Tin Electrode	99000087
Ag/AgCl Electrode	99000088
Gold Electrode	99000089
Ag/AgCl Sintered Electrode	99000058

DISPOSABLE ELECTRODES - EEG, EOG, ECG



Ambu™ Center Snap ECG Electrodes

Disposable ECG center snap electrode, white foam type, teardrop shaped with strong medical grade adhesive. Each package contains 5 disposable electrodes.

Product Description	Part Number
Disposable ECG electrode 1.7" x 1.5"	00103434
Disposable ECG electrode 1.9" x 2"	00103435



Softy Trode Repositionable ECG Electrode

Disposable ECG center snap electrode, cloth rectangular type. Electrode measures 1 1/4" x 7/8" and is packaged 6 electrodes per package, 450 electrodes per case.
Sold by the case.

Product Description	Part Number
Disposable Cloth ECG Electrode	00103337



Conmed Center Snap Cleartrace™ Electrodes

Disposable ECG Center Snap Electrode, The original clear tape electrode with conductive adhesive gel.

Product Description	Part Number
Disposable ECG - 5 pack	00101525
Cleartrace Electrode - 30 pack	00101526



Meditrace™ Mini Electrode

1 1/4" Teardrop shaped foam pre-gelled electrode. Smaller size makes this electrode ideal for facial or pediatric applications. Electrodes are packaged 30 per bag and have 20 bags per case. Sold per bag.


Product Description	Part Number
Medi-Trace Mini Electrode	5400-0012-00



Ambu™ Neuroline 720 Wet Gel Electrodes

Ambu Neuroline 720 electrodes are pre-gelled and self adhesive. They are ideal for use in sleep studies because of their small size and skin friendly adhesive that ensures the electrode stays in place for the entire examination even if the patient sweats. Each electrode comes with pre-attached leadwire that terminates into a 1.5mm connector.

Product Description	Part Number
Ambu 720 with 150cm with (60") leadwire - 10 pack	00103427
Ambu 720 Snap Electrode - 25 pack.	00103425



Ambu™ Single Patient Ground Electrode

Features unique "duck foot" design so that the corners cannot lift off. The concave sides ensure the electrode conforms to the contours of the body. The special green tab makes it easy to grasp the electrode without touching the adhesive area. Sold individually. 1.5mm connector.

Product Description	Part Number
Ambu Disposable ground with 100cm (40") Leadwire	00103465
Ambu Disposable ground with 150cm (60") Leadwire	00103466



Ambu™ Neuroline Solid Gel Electrodes

Ambu solid gel surface electrodes are self-adhesive and are both easy to reposition and comfortable for patients. The solid gel, combined with an effective adhesive allow the technician to locate the best recording site without leaving messy gel deposits. The special "duck-foot" design helps keep electrodes securely in place. Great for nerve conduction studies. Each package contains 12 electrodes per pouch. Ambu 700 solid gel electrodes have a 20mm x 15mm recording area.

Product Description	Part Number
Ambu 700 with 50cm (20") lead wire	00103433
Ambu 700 with 100cm (40") lead wire	00103432
Ambu 700 with 150cm (60") lead wire	00103431

CABLES - ECG LEADWIRES, EXTENSION JUMPERS, ADAPTORS



Snap Electrode Wires, 1.5mm connector

For use with snap type electrodes. Multicolor wire for easy identification. Super flexible and lightweight.

Each Package contains 10 leadwires.

Product Description	Part Number
10 Pack 48" Leadwire	00103468
10 Pack 60" Leadwire	7003-0001-00
10 Pack 72" Leadwire	00103469
10 Pack 96" Leadwire	00103470
5 Pack 26" Leadwire	7000-0198-01
2 Pack 72" Leadwire	7002-0004-00



Nonin Extension Cables

Compatible with Nonin reusable and disposable pulse oximeter probes. Sold individually.

Product Description	Part Number
1 Meter Extension Cable.....	00103473
3 Meter Extension Cable	00103474



Leadwire Extension Cable

Easily extend any cable using 1.5mm connectors with this cable.

Product Description	Part Number
48" Extension Cable	7002-0021-00



Jumper Cables

Stackable, 1.5mm connector. Sold in packs of 3

Product Description	Part Number
4" Stackable Jumpers.....	7002-0003-00
8" Stackable Jumpers	7002-0003-01



Universal 2mm and 1.5mm Adaptors

Adapt your 2mm pin to 1.5mm socket or 2mm socket to 1.5mm socket. Sold individually.

Product Description	Part Number
2mm pin to 1.5mm socket	97000121
2mm socket to 1.5mm socket.....	97000122