

PRESSURE TRANSDUCERS, CANNULA AND ACCESSORIES



■ Pro-Tech® PTAF2 Pressure Transducer

Can be used in either AC or DC configurations and during CPAP titrations. Each unit includes a 1075 Airflow cable, 1075 Snore Cable, 9V battery, Cannula Sample Pack, and PTAF2 sensor module. Please specify what type of recorder you are using to receive the correct interface cables.

Product Description

Pro-Tech PTAF2 Starter Kit (1291)00103584

Part Number



■ Pro-Tech® PTAF lite

Small, lightweight AC device, used during diagnostic studies. Dual output airflow and snore signal. Kit contains Cannula sample pack and PTAF lite device. 1.5mm Connector.

Product Description

Pro-Tech PTAF lite Starter Kit (1293) 00103588

Part Number



■ Salter Labs® BiNAPS™ Pressure/Snore Transducer

Detects airflow pressure and snore changes. Both airflow and snore channels feature Hi-Lo selectors to optimize signal levels. 1.5mm C

Product Description

Compumedics Pressure/Flow Module Kit – Siesta/Safiro 9007-0012-01

Part Number



■ Compumedics® AC Pressure/Flow Module

Detects respiratory airflow or CPAP pressure during titrations. Only compatible with Compumedics systems.

Product Description

BiNAPS Pressure/Snore Transducer (5500) 00103592

Part Number




■ Braebon ValueLine™ Pressure Sensor

Two output signals: Airflow and snoring. Requires no external power supply and is compatible with most recording systems. 1.5mm Connector.

Product Description

Braebon Pressure Sensor (5581) 00103589

Part Number



■ Pro-Tech® Pro-Flow™ and Pro-Flow Plus™ Cannulas

Pro-Tech ProFlow cannulas are designed for pressure transducers and provide optimal performance and comfort. They have soft, narrow prongs for minimal nasal resistance. Each cannula has a hydroponic, antimicrobial filter. Compatible with most pressure transducers.

Pro-Flow Nasal Cannula

Product Description

Adult Pro-Flow Nasal Cannula - 60 per case (1259) 00103590
Pediatric Pro-Flow Nasal Cannula - 60 per case (1263) 00103591

Part Number

Pro-Flow Plus Nasal Oral Cannula

Product Description

Adult Nasal/Oral Cannulas - 30 per case (1295) 5400-0020-00
Pediatric Nasal/Oral Cannulas - 30 per case (1354) 00103594

Part Number




Braebon Pressure Transducer Cannulas

Latex free nasal and nasal/oral cannula. Compatible with most pressure transducers. Each Cannula includes a safety filter.

Nasal Cannula

Product Description	Part Number
Value Flow Adult micro nasal cannula - 40 per case (5582)	00103595



Salter Labs® Sleep Diagnostic Cannula

Salter-Style ultra lightweight cannulas are designed to deliver clear signals for real time pressure/snore monitoring. Salter Sleep Diagnostic cannulas are compatible with most pressure transducers.

Nasal Cannula

Product Description	Part Number
Adult Pressure/Snore Cannula with hydrophobic/anti-microbial filter Male Luer Lok® connector 25/case (5011)	00103601
Adult Pressure/Snore Cannula w/o filter Male/Female Luer Lok® connector 25/case (5012)	00103688
Pediatric Pressure/Snore Cannula with hydrophobic/anti-microbial filter Male Luer Lok® Connector 25/case (5021)	00103689
Pediatric Pressure/Snore Cannula w/o filter Male/Female Luer Lok® connector 25/case (5022)	00103690



Nasal/Oral Cannula

Product Description	Part Number
Adult Nasal/Oral Airflow Pressure/Snore Cannula with non-ferrous form wire in oral prong and hydrophobic/anti-microbial filter Male Luer Lok® connector 25 per case (5001)	00103691
Adult Nasal/Oral Airflow Pressure/Snore Cannula with non-ferrous form wire in oral prong w/o hydrophobic/anti-microbial filter Male/Female Luer Lok® connector 25 per case (5002)	00103692



Nasal ETCO2 Cannula

Product Description	Part Number
Adult Nasal CO2 Sample Cannula with 7' line Female Luer Lok® connector 25 per case (4000F)	5400-0011-01
Adult Nasal CO2 Sample Cannula with 2' line Female Luer Lok® connector 25 per case (4000F-2)	5400-0021-02
Pediatric Nasal CO2 Sample Cannula with 7' line Female Luer Lok® connector 25 per case (4100F)	00103693



Braebon Pressure Transducer Cannulas

Latex free nasal and nasal/oral cannula. Compatible with most pressure transducers. Each Cannula includes a safety filter.

Product Description	Part Number
For oxygen flow up to 6 LPM Adult Clear Cannula with 7' supply tube 50 per case	00103553
For oxygen flow up to 3 LPM Adult "Micro" Cannula with 7' supply tube 50 per case	00103694



Braebon Pressure Transducer Cannulas

Latex free nasal and nasal/oral cannula. Compatible with most pressure transducers. Each Cannula includes a safety filter.

Product Description	Part Number
Male to Male Luer Adaptor	
- Sold individually	00103356
Female to Female Adaptor	
- 10 Pack	7000-0145-00