

## PULSE OXIMETRY



### ■ Nonin® Pulse Oximeter

Available in both 3011 and 3012 models for use with Compumedics E-Series, Siesta, and Somte sleep diagnostic equipment. Complete kit includes; oximeter, 3m probe, and one package of adult flexi-wrap. The 3012 model provides the additional output of the plethysmograph waveform.

Product Description	Part Number
3011 Complete Kit .....	7000-0238-01
3012 Complete Kit .....	7000-0239-01
3011 Oximeter ONLY .....	7008-0030-01
3012 Oximeter ONLY .....	7008-0040-01



### ■ Nonin® Disposable Oximeter Sensors

Compatible with Nonin pulse oximeters. 10 sensors per pack.

Product Description	Part Number
Adult Flexi-form sensor w/40" leadwire .....	00103463
Pediatric Flexi-form sensor w/40" leadwire .....	00103464




7001-0002-00

### ■ Nonin® Pulse Oximeter Sensors

Replacement oximeter probes and accessories for Nonin oximeters.

Product Description	Part Number
Reusable 1M Finger Probe (8000J) includes flexi-wrap .....	7001-0002-00
Reusable 3M Finer Probe (8000J-3) includes flexi-wrap .....	7001-0008-00
Adult Flex-wrap 25 per pack .....	7000-0147-01
Reusable 1 Meter Adult Articulated Finger Clip .....	00103716
Reusable 3 Meter Adult Articulated Finger Clip .....	00103717
Reusable 1 Meter Pediatric Articulated Finger Clip .....	00103718
Reusable 3 Meter Pediatric Articulated Finger Clip .....	00103719
3 Foot Extension Cable .....	00103473
9 Foot Extension Cable .....	00103474
10 Foot Extension Cable – P-Series .	3101-0003-00



7000-0147-01



00103716

### ■ Nonin Soft Sensors

Compatible with Nonin pulse oximeters. Available in 3 convenient sizes. Latex-Free.

Product Description	Part Number
1 Meter Small - 7.5mm to 12mm .....	00103740
1 Meter Medium - 10mm to 19mm .....	00103741
1 Meter Large - 12.5mm to 25.5mm .....	00103742
3 Meter Small - 7.5mm to 12mm .....	00103747
3 Meter Medium - 10mm to 19mm .....	00103748
3 Meter Large - 12.5mm to 25.5mm .....	00103749

