

RESPONSEPixx™

Ergonomic response boxes



FEATURES

- ERGONOMIC DESIGN
- ILLUMINATED BUTTONS CAN CUE SUBJECT RESPONSES
- TTL OUTPUTS ELIMINATE LATENCY INHERENT IN USB SOLUTIONS
- MICROSECOND-PRECISE RESPONSE TIMES WHEN USED WITH DATAPixx OR VIEWPixx
- RELIABLE SOLUTION RECORDS ALL BUTTON PRESSES
- CUSTOM LAYOUTS AVAILABLE

OVERVIEW

There are 3 versions of our RESPONSEPixx line of reliable, ergonomic response boxes. Each of the buttons can be illuminated under software control, and can cue subjects during forced-choice paradigms. The RESPONSEPixx uses medical grade silicone cables, resulting in maximum flexibility and durability. Plugging the RESPONSEPixx into a DATAPixx or VIEWPixx will give you microsecond-precise response times, with a reliability that exceeds USB-based response boxes. Custom configurations are available.

RESPONSEPixx SPECIFICATIONS

RESPONSEPixx HANDHELD

P/N: VPX-ACC-3100

SPECIFICATIONS:

- Handheld response box with 5 thumb buttons
- Guaranteed microsecond-precise reaction time measurement when interfaced to DATAPixx system
- Illuminated buttons can cue subjects when a response is required, and which responses are valid

DIMENSIONS:

14 x 10.2 x 3.5 cm (5.5 x 4.0 x 1.38 inches)
3.66m (12 ft) interface cable with DB-25 connector



RESPONSEPixx DUAL HANDHELD

P/N: VPX-ACC-3000

SPECIFICATIONS:

- Two handheld response boxes with 2 thumb buttons each
- Guaranteed microsecond-precise reaction time measurement when interfaced to DATAPixx system
- Illuminated buttons can cue subjects when a response is required, and which responses are valid

DIMENSIONS:

10 x 5.08 x 2.85 cm (3.95 x 2.0 x 1.12 inches)
3.66 m (12 ft) interface cable with DB-25 connector



RESPONSEPixx TABLETOP

P/N: VPX-ACC-3050

SPECIFICATIONS:

- Tabletop response box with 5 large robust buttons
- Guaranteed microsecond-precise reaction time measurement when interfaced to DATAPixx system
- Illuminated buttons can cue subjects when a response is required, and which responses are valid

DIMENSIONS:

23.3 x 13.34 x 6 cm (9.18 x 5.25 x 2.36 inches)
3.66 m (12 ft) interface cable with DB-25 connector



VPixx Technologies Inc.

1494 Montarville suite 206
Saint-Bruno, QC
Canada, J3V 3T5

TEL/FAX: (514) 328-7499
EMAIL: sales@vpixx.com

